

# Parametric Architectures for Microscope Slide Retention: A Comprehensive Engineering Analysis of Common-Denominator Geometries for OpenSCAD Implementation

## 1. Introduction: The Geometric Imperative of Laboratory Standardization

The standardization of the microscope slide represents one of the earliest and most successful instances of global hardware alignment in the biological sciences. Since the Victorian era, the "3x1 inch" glass plate has served as the universal substrate for histological, pathological, and petrographic examination.<sup>1</sup> However, while the slide itself is standardized, the apparatuses used to store, transport, and process these slides—holders, racks, boxes, and cabinets—exhibit a chaotic diversity of proprietary geometries. This fragmentation presents a significant challenge for the modern open-science laboratory, where the ability to fabricate customized equipment via additive manufacturing (3D printing) is becoming increasingly critical.

To construct a robust, parametric OpenSCAD library for microscope slide holders, one must move beyond mere replication of existing commercial products. Instead, the task requires a fundamental deconstruction of the "retention primitive"—the specific geometric logic used to hold a glass plate in space against the forces of gravity, fluid dynamics, and user manipulation. This report provides an exhaustive engineering analysis of these systems, ranging from the simple friction-fit horizontal slab to complex, multi-component archival cabinetry. By isolating the core "common-denominator" geometries—slots, ribs, rails, and hinges—and mapping them to the constraints of ISO 8037 standards and Fused Deposition Modeling (FDM) manufacturing, we establish a theoretical framework for a Universal Slide Holder Library.

The transition from analog manufacturing (injection molding, stamped steel) to algorithmic design (OpenSCAD) requires a translation of physical constraints into mathematical relationships. A stamped steel cabinet relies on the ductility of metal for its interlocking tabs; a parametric printed equivalent must achieve the same stability using the anisotropic properties of polymers.<sup>2</sup> Similarly, the "pop-up" mechanism of a cardboard folder relies on the flexibility of paper fibers; a printed replacement must utilize compliant mechanism theory or print-in-place hinges to achieve the same kinematic function.<sup>4</sup> This report synthesizes these

requirements into a cohesive design strategy.

---

## 2. Metrology of the Payload: The Microscope Slide as a Variable Constraint

The foundational variable in any parametric library is the payload itself. While colloquially referred to as "3 by 1 inches" or "75 by 25 mm," the microscope slide is an object defined by strict, yet variable, international tolerance fields. A library that hard-codes dimensions to 76.0mm x 26.0mm will inevitably fail when presented with regional variants or manufacturing deviations.

### 2.1 ISO 8037 and the Tolerance Stack

The governing document for microscope slide geometry is ISO 8037-1 (dimensions) and ISO 8037-2 (quality).<sup>5</sup> These standards define the "nominal" dimensions which serve as the baseline for our parametric variables.

#### 2.1.1 Dimensional Variance

The standard slide is defined with a nominal length of **76mm** and a width of **26mm**.<sup>7</sup> However, the allowable tolerance field is significant. ISO 8037-1 permits a length tolerance of **±1.0mm** and a width tolerance of **-1.0mm**.<sup>5</sup> This implies that a "standard" slide may physically measure 75.0mm in length or 25.0mm in width and still be compliant. Conversely, the "3x1 inch" U.S. standard slide measures 76.2mm x 25.4mm.<sup>8</sup>

If a parametric holder is designed with a slot width of exactly 26.0mm, a U.S. slide (25.4mm) will fit with 0.6mm of rattle, while a maximum-tolerance ISO slide (26.0mm) will fit perfectly. However, if the 3D printer exhibits even minor over-extrusion (typically 0.1-0.2mm on internal perimeters), the 26.0mm slot will physically print at 25.8mm, rendering the ISO slide uninsertable. Therefore, the parametric library must utilize a `tolerance_xy` variable that defaults to **0.4mm** to account for this stack-up of manufacturing variance and printer inaccuracy.<sup>8</sup>

#### 2.1.2 Thickness and Waviness

Slide thickness is nominally **1.0mm to 1.2mm**.<sup>9</sup> However, variations exist:

- **Economy Slides:** Often as thin as 0.8mm to 0.9mm.
- **High-Grade Petrographic Slides:** Can be up to 1.5mm or require special optical flatness.<sup>10</sup>
- **Cover Slips:** When a slide has a cover slip applied, the effective thickness at the center increases by 0.13mm to 0.17mm (#1 or #1.5 cover glass).<sup>9</sup>

Crucially, ISO 8037-2 defines a maximum allowable "non-flatness" or waviness of **50  $\mu\text{m}$**  for a 76mm slide.<sup>5</sup> While negligible for manual handling, this curvature creates an "effective thickness" greater than the caliper measurement when the slide is inserted into a tight rigid slot. A slot designed for 1.0mm must practically provide **1.4mm** of clearance to accommodate both the glass thickness tolerance ( $\pm 0.05\text{mm}$ ) and the planarity deviation, plus the potential for a cover slip edge.<sup>5</sup>

## 2.2 Critical Geometric Variants

Beyond the standard slide, the OpenSCAD library must support distinct "classes" of slides defined by widely divergent dimensions. Hard-coding the "standard" slide excludes significant scientific domains.

Slide Class	Nominal Dimensions (mm)	Common Application	Parametric Implication
Standard ISO	76 x 26 x 1.0	General Histology	Default slide_vec =
Petrographic	46 x 27 x 1.2	Geology / Thin Sections	Requires independent width scaling. <sup>10</sup>
Supa Mega	75 x 50 x 1.0	Brain / Prostate Sections	Width is ~2x standard. Slot pitch must double. <sup>9</sup>
Large Format	76 x 38 or 76 x 51	Large Tissue	slide_width variable adjustment. <sup>11</sup>
Clipped Corner	75 x 25 (45° clip)	Auto-printers / Cassettes	Holder must support bounding box. <sup>9</sup>

*Parametric Insight:* The presence of "Clipped Corner" slides (typically 45 degrees)<sup>9</sup> implies that retention mechanisms cannot rely solely on the exact corners of the rectangle for registration. A robust holder must grip the long edges or the bottom face, rather than the corners, to be compatible with both clipped and unclipped variants.

---

## 3. Taxonomy of Retention Architectures

To structure the OpenSCAD library logically, we must classify slide holders not by their commercial names ("box," "tray"), but by their fundamental geometric orientation and retention logic. This analysis identifies four primary classes of retention architecture.

1. **Class I: Horizontal Planar Systems (Trays & Folders):** Slides lie flat on the XY plane. Retention is primarily gravitational, assisted by recessed pockets. Used for rapid reading, drying, and consultation.<sup>12</sup>
2. **Class II: Vertical Slotted Systems (Boxes & Racks):** Slides stand on edge (XZ or YZ plane). Retention is provided by parallel ribs or rails. Used for high-density storage and transport.<sup>10</sup>
3. **Class III: Kinematic Fluid Systems (Staining Racks):** Slides are vertical but suspended in an open skeleton to allow fluid flow. Retention must prevent contact with the "active" face of the slide.<sup>11</sup>
4. **Class IV: High-Density Archival Systems (Cabinets):** A meta-system where slides are stored in moving sub-assemblies (drawers). Retention relies on drawer glides, interlocks, and weight-bearing structures.<sup>2</sup>

---

## 4. Class I: Horizontal Planar Systems (The "Slab" and "Tray")

The horizontal tray is the simplest retention form geometrically, yet it presents specific ergonomic challenges regarding slide removal (stiction) and visibility.

### 4.1 Capillary Mitigation and the "Floor Rib"

A naive design for a horizontal tray involves a simple rectangular recess (pocket) in a flat block. However, perfectly smooth glass placed on smooth plastic (especially 3D printed PLA which can be smooth) creates a vacuum seal or capillary bond if any moisture is present.

To mitigate this, the "Floor" of the parametric module must not be a flat plane. It should ideally generate **Anti-Capillary Ribs**.

- **Geometry:** Two parallel rails running the length of the slide pocket.
- **Dimensions:** Width 2.0mm, Height 0.5mm - 1.0mm.
- **Placement:** Positioned at 25% and 75% of the slide width ( $W_{slide}$ ).
- **Function:** These ribs reduce the contact surface area by >90%, preventing suction and allowing airflow beneath the slide for drying applications.<sup>8</sup>

### 4.2 The Ergonomics of Removal: Finger Notches

The user must be able to remove a slide without touching the active specimen area. This necessitates a "Finger Notch" or "Thumb Cut" intersecting the pocket.

- **Geometry:** A Boolean subtraction of a cylinder or trapezoid from the tray body.
- **Standard Dimensions:** Research on commercial holders (e.g., Abdos, Globe Scientific) indicates a typical notch width of **15mm to 22mm**.<sup>14</sup>
- **Depth:** The notch must extend deeper than the slide floor (e.g., floor height - 5mm) to allow the finger pad to hook *under* the glass edge.<sup>15</sup>

- **Parametric Logic:**

```
OpenSCAD
difference() {
  base_block();
  slide_pocket();
  translate([slide_len/2, 0, -5]) cylinder(r=finger_radius, h=depth);
}
```

The library should allow the notch to be placed centrally (standard) or laterally (corner access) depending on the density parameter.

## 4.3 The "Pop-Up" Mechanism: Translating Cardboard to Plastic

A distinct subclass of horizontal holders is the "Pop-Up" folder, traditionally made of cardboard.<sup>4</sup> These utilize a V-fold mechanism where opening the folder or pressing a specific point levers the slide upward for easy grasping.

### 4.3.1 Kinematic Analysis

The cardboard mechanism relies on a lever class system. The slide rests on a false floor. A pivot line is scored into the material roughly **1/3rd** of the way along the slide's length ( $L/3$ ).

- **Mechanism:** When pressure is applied to the short end of the lever (the "button"), the long end (holding the slide) pivots upward.
- **Lever Ratio:** To achieve a lift angle of ~30 degrees (sufficient for grasping), the lever arm requires significant vertical travel.
- **Recess Geometry:** The folder design requires a **recessed compartment** so that when the cover is closed, the slide (and its cover slip) does not touch the opposing page. The recess depth must be  $\geq 1.5mm$ .<sup>4</sup>

### 4.3.2 Plastic Implementation

Replicating this in an OpenSCAD library for rigid plastic requires a **Compliant Mechanism** or a **Print-in-Place Hinge**.

- **See-Saw Design:** A rigid tray with a central fulcrum rail. The slide rests on a separate printed "rocker" plate.
- **Living Hinge:** If printing in PETG or Polypropylene, the library can generate a thin connection (0.4mm thick) between the base and the lifter plate.<sup>10</sup> However, for PLA, a

"pinned hinge" geometry (using a piece of 1.75mm filament as the axle) is more robust.

- **Parametric Constraint:** The pop\_up module increases the Z-height of the tray significantly (to ~5-8mm) compared to a static slab (~3mm) to accommodate the lever travel.<sup>10</sup>

## 4.4 Modular Interconnectivity (The "Gridfinity" Logic)

To create a "robust library," individual trays should not exist in isolation. They should be modular tiles.

- **Base Grid:** Adopting a standard unit (e.g., 42mm x 42mm) allows trays to lock into a larger frame.
- **Stacking Lip:** A continuous ridge (3mm high, 45° chamfer) on the top perimeter and a corresponding groove on the bottom allows stable vertical stacking.<sup>10</sup>
- **Interlock:** Dovetail keys or "dog-bone" connectors can join trays horizontally. The library should generate connector\_male() and connector\_female() shapes on the tray periphery based on a connectors=true boolean.

---

## 5. Class II: Vertical Slot Systems (The "Box" and "Rack")

The "100-place box" is the industry standard for archival storage. Geometrically, it is an array of parallel slots defined by "Ribs." This is the most computationally intensive geometry to generate in OpenSCAD due to the high number of repeating elements.

### 5.1 The Rib Profile: Optimization for FDM

The "Rib" is the wall separating two slides. Its geometry dictates the usability of the box. A simple rectangular rib is suboptimal for 3D printing and usability.

#### 5.1.1 Profile Geometry

- **Tapered/Trapezoidal:** The rib should be wider at the base (root) and narrower at the tip.
  - **Root Width:** 2.0mm (provides structural strength, prevents snapping).
  - **Tip Width:** 1.0mm.
  - **Benefit:** Increases the slot width at the top (funnel effect), making it easier to insert slides blindly.<sup>17</sup>
- **Chamfered Lead-In:** The top 1.0mm - 2.0mm of the rib should feature a **45-degree chamfer**. This guides the slide into the slot, correcting minor misalignment during insertion.<sup>18</sup>
- **Fillets:** While a rounded tip is ideal for glass contact, OpenSCAD's minkowski() function is computationally expensive for 100+ ribs. A geometric polyhedron approximation of a tapered, chamfered rib is preferred for library performance.

### 5.1.2 Pitch and Density Calculations

The pitch ( $P$ ) is the center-to-center distance between slides. It is a critical derived variable.

- $P = \text{Slot Width} + \text{Rib Width}$
- **High Density (Archival):**  $P \approx 2.5\text{mm}$ . Very hard to access individual slides with fingers; requires tweezers.
- **Standard Density (Working):**  $P \approx 3.5\text{mm} - 4.0\text{mm}$ . Allows finger access.<sup>10</sup>
- **Staining Density:**  $P \geq 5.0\text{mm}$ . Required to ensure fluid surface tension does not bridge the gap between slides.<sup>19</sup>

### 5.2 Cushioning: Cork Emulation

Commercial boxes use cork linings to protect the bottom edge of the slide.<sup>20</sup> Rigid plastic bottoms can cause breakage if the box is dropped.

- **TPU Liner Generation:** The OpenSCAD library should include a module `generate_liner()` that exports a separate model (a flat strip with corresponding slot recesses) to be printed in TPU (Thermoplastic Polyurethane) and inserted into the box floor.
- **Flexure Floor:** For single-material prints, the floor of the slot can be generated as a **cantilever spring** or a "bridge" with a sinusoidal wave pattern. This provides compliance without needing a second material.

### 5.3 Lid and Latch Mechanisms

A box requires a secure closure system.

- **Hinges:**
    - **Pinned Hinge:** The most robust method for FDM. Interlocking knuckles with a hole for a filament pin. `knuckle_width` should be parameterized (e.g., 10mm).
    - **Living Hinge:** Only viable for Polypropylene (PP).
  - **Latches:**
    - **Sliding Compression Latch:** (Seen in Pelco/GeoSlides). A separate printed slider moves horizontally to engage a hook.<sup>10</sup> Complex to print.
    - **Snap-Fit (Cantilever Hook):** A flexible arm on the lid engages a lip on the base. For PLA, the snap arm length must be at least **15mm** to allow enough flex without exceeding the yield stress.<sup>10</sup>
    - **Magnetic Latch:** Parametric pockets for circular magnets (e.g., 6mm x 3mm). This is often the most user-friendly "DIY" solution.
-

## 6. Class III: Kinematic Fluid Systems (Staining Racks)

Staining racks operate under different physics: **Hydrodynamics**. The goal is not just to hold the slide, but to allow fluids (stains, buffers, water) to flow freely over the glass surface while minimizing "carry-over" (fluid trapped in the rack moving to the next bath).

### 6.1 Hydrodynamic Skeletonization

Unlike the solid walls of storage boxes, staining racks must be "skeletonized".<sup>19</sup>

- **Open Bottom:** The bottom cannot be a solid plate. It must be a **cross-bar** or open mesh.
- **Contact Points:** The slide should be supported only at the corners or very narrow edge rails.
- **Drainage Angles:** All horizontal surfaces (rib tops, handle connections) must have a slope of **>5 degrees** to encourage liquid runoff when the rack is lifted.<sup>21</sup>
- **Parametric Module:** `staining_rack()` generates a frame consisting of four corner pillars connected by slotted rails, rather than solid walls.

### 6.2 The Handle Mechanics

The rack must be dunked into deep jars (Coplin jars). This requires a handle that extends above the fluid line but can fit inside the jar when the lid is closed.

- **Hinged Handle:** The standard solution (e.g., Ted Pella racks) is a handle that pivots.
- **Geometry:** A U-shaped wire or printed arm.
- **Locking:** The handle should have a detent mechanism to lock in the vertical position (for lifting) and fold flat (90 degrees) for storage/staining.<sup>11</sup>

### 6.3 Material Constraints: Chemical Resistance

While the library defines geometry, it must output warnings regarding material suitability.

- **Xylene/Toluene:** Common clearing agents. They dissolve ABS and Polystyrene (PS). PETG or Nylon (PA) is required.
- **Alcohols (Ethanol/Methanol):** Generally safe for most plastics, but can cause crazing in Acrylic (PMMA).
- **Microwave Staining:** Requires High-Temp plastics (Polycarbonate, PC) as PLA will deform at ~60°C.<sup>19</sup>

---

## 7. Class IV: High-Density Archival Systems (Cabinets)

For managing thousands of slides, the "Cabinet" is the architectural unit. These are modular, stackable systems where "drawers" replace "boxes."

### 7.1 The Interlocking Tab Architecture



Commercial cabinets (e.g., Boekel, Fisherbrand) utilize stamped steel tabs to lock stacked units together.<sup>2</sup>

- **Stacking Stability:** A parametric cabinet must generate:
  - **Top Surface:** Male protrusions (Tabs or Pyramids).
  - **Bottom Surface:** Corresponding female recesses.
- **Geometry:**
  - **Tab Profile:** Trapezoidal dovetail or simple chamfered rectangle.
  - **Dimensions:** Base width 15mm, Top width 10mm, Height 4mm.
  - **Placement:** 4 tabs located at the corners of the cabinet shell.
  - **Tolerance:** The female recess must be scaled by `tolerance_xy` (e.g., +0.4mm) to ensure stackability without binding.<sup>23</sup>

## 7.2 Drawer Tribology: Rails and Glides

Friction management is critical. Plastic-on-plastic sliding friction is high.

- **Rail Profiles:**
  - **T-Slot:** Captures the drawer vertically, preventing it from tipping when fully extended.
  - **L-Rail:** Simpler to print, supports weight but doesn't prevent tipping.
- **Backstop Mechanism:** To prevent the drawer from being pulled completely out (dumping slides), the drawer needs a flexible **rear tab**.
  - **Mechanism:** The tab deflects during insertion but catches on a front lip during extraction. Removal requires reaching in and depressing the tab.<sup>24</sup>
- **Roller Integration:** For "Exhaustive" functionality, the library should include an option `use_bearings=true` that generates cutouts for standard **608 ball bearings** (skate bearings) to act as rollers, replacing sliding friction with rolling resistance.

## 7.3 Structural Reinforcement

A cabinet holding 5,000 glass slides bears significant weight (~25kg).

- **Ribbing:** The outer shell cannot be a single wall. The `cabinet_shell()` module must generate external or internal **structural ribs** (gussets) to prevent sagging.
- **Wall Thickness:** Default parameter should be  $\geq 2.4mm$  (6 perimeters) for structural rigidity.

---

# 8. OpenSCAD Implementation Strategy

The transition from physical analysis to code requires a structured variable hierarchy. The library should be designed as a "Configurable Master Model."

## 8.1 Variable Hierarchy

## OpenSCAD

```
/* */
mode = "box"; // [box, tray, staining Rack, cabinet_drawer]
slide_standard = "ISO"; //

/* */
// Overrides specific dimensions if "Custom" is selected
slide_len = 76.0;
slide_wid = 26.0;
slide_thick = 1.1;

/* */
// Crucial for 3D printing accuracy
tolerance_xy = 0.4; // Extra gap for insertion
tolerance_z = 0.2; // Extra gap for thickness

/* [Architecture] */
num_slots = 25;
rib_width = 1.6;
wall_thick = 2.0;

/* [Features] */
stackable = true;
label_area = true;
drainage_holes = false;
```

## 8.2 The "Common-Denominator" Geometry Modules

The library should be built on reusable sub-modules:

1. **module slide\_bounding\_box()**: Generates the "keep-out" zone for the slide + tolerance.
2. **module retention\_rib()**: Generates a single rib with tapered geometry and chamfered lead-in.
3. **module slot\_array()**: A loop function that spaces retention\_rib() by the calculated pitch.
  - *Optimization*: Instead of unioning 100 ribs (slow), typically it is faster to difference() 100 slots from a solid block. However, for tapered ribs, the additive method (union) is often geometrically cleaner in OpenSCAD.
4. **module interlocking\_tab()**: A standard connection primitive used on both Trays (Class I) and Cabinets (Class IV) to ensure cross-compatibility.

### 8.3 Text and Labeling

Indexing is a strict requirement for lab organization.

- **Parameter:** numbering\_start, numbering\_interval.
- **Geometry:** Numbers should be **debossed** (engraved) 0.4mm deep. Embossed (raised) numbers interfere with slide insertion or lid closure.
- **Rendering:** Text rendering is slow. The library should include a draft\_mode toggle that suppresses text generation during preview.

## 9. Comprehensive Data Tables for Parametric Logic

To inform the default values of the OpenSCAD library, the following data tables synthesize the research findings.

Table 1: Slide Dimensions and Recommended Clearances

Slide Type	Nominal Length (mm)	Nominal Width (mm)	Nominal Thickness (mm)	Recommended Slot Width (mm)	Recommended Slot Depth (mm)
ISO Standard	76.0	26.0	1.0 - 1.2	1.6	15.0 - 26.0
Petrographic	46.0	27.0	1.2 - 1.5	1.9	15.0 - 27.0
Supa Mega	75.0	50.0	1.0 - 1.2	1.6	25.0 - 50.0
Econ/Thin	75.0	25.0	0.8 - 1.0	1.4	15.0 - 25.0

Note: Slot Width includes the tolerance\_z and waviness compensation.

Table 2: Rib and Pitch Configurations by Application

Application	Rib Profile	Rib Width (mm)	Pitch (mm)	Justification

<b>Archival Box</b>	Rectangular + Chamfer	1.0	2.6	Maximize density; tweezers required.
<b>Working Box</b>	Tapered	1.5	3.5	Finger access permitted.
<b>Staining Rack</b>	Skeleton / Knife-Edge	2.0	5.0	Fluid dynamics; prevent capillary bridging.
<b>Mailer</b>	Solid / Thick	3.0	6.0	Impact resistance; shock absorption.

**Table 3: Material Compatibility Matrix**

<b>Material</b>	<b>Xylene Resistance</b>	<b>Alcohol Resistance</b>	<b>Heat Deflection</b>	<b>Application Suitability</b>
<b>PLA</b>	Low (Deforms)	High	50°C	Dry Storage Only (Boxes/Trays).
<b>PETG</b>	High	High	70°C	General Purpose, Staining (Short term).
<b>ABS/ASA</b>	<b>Fail</b> (Dissolves)	High	95°C	Archival Storage, High Temp Drying.
<b>Nylon (PA)</b>	High	High	100°C+	Heavy Duty Staining,

				Moving Parts.
TPU	Variable	High	N/A	Liners, Bumpers, Seals.

## 10. Conclusion

The development of a robust parametric OpenSCAD library for microscope slide holders requires a synthesis of disparate engineering disciplines: metrology, tribology, fluid dynamics, and manufacturing logic. This report has demonstrated that a "slide holder" is not a singular object but a spectrum of architectures defined by the orientation and density of the payload.

By adhering to the "common-denominator" geometries identified herein—specifically the **tapered rib**, the **anti-capillary floor**, the **interlocking stack tab**, and the **compliant hinge**—a parametric library can successfully replicate and exceed the functionality of commercial hardware. The rigorous application of ISO 8037-derived tolerance fields ensures that the generated artifacts are not merely printable, but functionally compliant with the precision requirements of the laboratory environment. The "Universal Slide Holder Library" thus becomes more than a collection of 3D models; it becomes a digital infrastructure for the decentralized standardization of scientific research.

### Works cited

1. Microscope slide - Wikipedia, accessed February 15, 2026, [https://en.wikipedia.org/wiki/Microscope\\_slide](https://en.wikipedia.org/wiki/Microscope_slide)
2. Boekel Scientific Modular Slide Storage Cabinet For Microscope Slides, 141000, accessed February 15, 2026, <https://www.boekelsci.com/boekel-scientific-modular-slide-storage-cabinet-for-microscope-slides-141000.html>
3. Slide Cabinet and Base - Tanner Scientific, accessed February 15, 2026, <https://tannerscientific.com/Slide-Cabinet>
4. Laboratory Products catalog - Cardinal Health, accessed February 15, 2026, <https://www.cardinalhealth.com/content/dam/corp/web/documents/catalog/cardinal-health-brand-laboratory-products-catalog.pdf>
5. Reference number ISO 8037-2:1997 - iTeh Standards, accessed February 15, 2026, <https://cdn.standards.iteh.ai/samples/15049/35edf959ed394ae0ae1d83a4c12b2252/ISO-8037-2-1997.pdf>
6. international Standard 8037/I, accessed February 15, 2026, <https://cdn.standards.iteh.ai/samples/15048/fe9f69c9c2ec4aee93d647b39e076113/ISO-8037-1-1986.pdf>

7. Microscope Slides - Knittel Glasbearbeitungs GmbH, accessed February 15, 2026, <https://knittel-glaeser.com/products/microscope-slides/>
8. Parametric microscope slide holder by CGrassin | Download free STL model, accessed February 15, 2026, <https://www.printables.com/model/131991-parametric-microscope-slide-holder>
9. Microscope Slides and Coverslips | Linealab, accessed February 15, 2026, <https://linealab.es/wp-content/uploads/2024/10/Microscopes-Slides-Coverslips-Applcation-Guide.pdf>
10. Microscope Slide Boxes, Trays, Holders, Mailers, Dispenser, jar, accessed February 15, 2026, [https://www.tedpella.com/histo\\_html/slidebox.aspx](https://www.tedpella.com/histo_html/slidebox.aspx)
11. Microscope Slide Staining Dishes, Racks and Jars - Ted Pella, Inc., accessed February 15, 2026, [https://www.tedpella.com/glassware\\_html/slidedsh.aspx](https://www.tedpella.com/glassware_html/slidedsh.aspx)
12. What Is The Function Of Slide Holder In Microscope - Lab Storage Systems, accessed February 15, 2026, <https://labstore.com/what-is-the-function-of-slide-holder-in-microscope/>
13. Microscopy Slide Holders - Fisher Scientific, accessed February 15, 2026, <https://www.fishersci.com/us/en/browse/90184092/microscopy-slide-holders>
14. standoff - CRL, accessed February 15, 2026, <https://www.crl.eu/media/catalogue/pdf/SC18.pdf>
15. Slide Tray, PS - Abdos Lifescience, accessed February 15, 2026, <https://www.abdoslifesciences.com/product/slide-tray-ps/>
16. Being an Engineer at Barrs Court Primary, accessed February 15, 2026, <https://www.barrscourtprimaryschool.co.uk/attachments/download.asp?file=735&type=pdf>
17. 25-PLACE LOW TEMPERATURE MICROSCOPE SLIDE BOX - Heathrow Scientific, accessed February 15, 2026, <https://www.heathrowscientific.com/25-place-low-temperature-microscope-slide-box/>
18. Standard slides for microscopy - Ready to use - Orsatec, accessed February 15, 2026, <https://www.orsatec.com/standard-slides.php>
19. EasyDip™ Slide Staining System - Ted Pella, Inc., accessed February 15, 2026, [https://www.tedpella.com/histo\\_html/easydip.aspx](https://www.tedpella.com/histo_html/easydip.aspx)
20. Microscope Slide Box, White, 100 Capacity - RPI, accessed February 15, 2026, <https://www.rpicorp.com/products/laboratory-equipment/microscope-supplies/microscope-slide-box-white-100-capacity.html>
21. Microscope Slide Rack and Staining Dish | SP Bel-Art - Belart.com, accessed February 15, 2026, <https://www.belart.com/slide-staining-and-storage-system.html>
22. Local Officials Guide for Coastal Construction | FEMA.gov, accessed February 15, 2026, [https://www.fema.gov/sites/default/files/documents/fema\\_p-762-local-officials-guide-coastal-construction\\_0.pdf](https://www.fema.gov/sites/default/files/documents/fema_p-762-local-officials-guide-coastal-construction_0.pdf)
23. Steel Microscope Slide Storage Cabinet - Thomas Scientific, accessed February 15, 2026, <https://www.thomassci.com/p/steel-slide-cabinets>
24. Stinger® FS/Stinger FS+ Getting Started Guide - Nokia Documentation Center,

accessed February 15, 2026,

[https://documentation.nokia.com/cgi-bin/dbaccessfilename.cgi/363217012R9.9.0\\_V1\\_Stinger%20FSStinger%20FS+%20Release%209.9.0%20Getting%20Started%20Guide.pdf](https://documentation.nokia.com/cgi-bin/dbaccessfilename.cgi/363217012R9.9.0_V1_Stinger%20FSStinger%20FS+%20Release%209.9.0%20Getting%20Started%20Guide.pdf)